Instruction Manual

HI 2214 HI 2215

pH/mV/°C Bench Meter with Custom Buffers and Interval Logging





Dear Customer,

Thank you for choosing a Hanna Instruments product.

Please read this instruction manual carefully before using these instruments. This manual will provide you with the necessary information for correct use of these instruments, as well as a precise idea of their versatility. If you need additional technical information, do not hesitate to e-mail us at tech@hannainst.com or view our worldwide contact list at www.hannainst.com.

WARRANTY

HI 2214 and **HI 2215** are guaranteed for two years against defects in workmanship and materials when used for their intended purpose and maintained according to instructions. Electrodes and probes are guaranteed for six months. This warranty is limited to repair or replacement free of charge.

Damage due to accidents, misuse, tampering or lack of prescribed maintenance is not covered.

If service is required, contact the dealer from whom you purchased the instrument. If under warranty, report the model number, date of purchase, serial number and the nature of the problem. If the repair is not covered by the warranty, you will be notified of the charges incurred. If the instrument is to be returned to Hanna Instruments, first obtain a Returned Goods Authorization number from the Technical Service department and then send it with shipping costs prepaid. When shipping any instrument, make sure it is properly packed for complete protection.

TABLE OF CONTENTS	
WARRANTY	2
PRELIMINARY EXAMINATION	3
GENERAL DESCRIPTION	3
FUNCTIONAL DESCRIPTION	
HI 2214 SPECIFICATIONS	
HI 2215 SPECIFICATIONS	
OPERATIONAL GUIDE	
pH CALIBRATION	9
PH BUFFER TEMPERATURE DEPENDENCE	
RELATIVE mV CALIBRATION	16
GOOD LABORATORY PRACTICE (GLP)	1/
LOGGING FUNCTION	
SETUP	26
TEMPERATURE CALIBRATION (for technical personnel only)	30
mV CALIBRATION (for technical personnel only)	
PC INTERFACE ELECTRODE CONDITIONING & MAINTENANCE	3Z
TROUBLESHOOTING GUIDE	
TEMPERATURE CORRELATION FOR pH SENSITIVE GLASS	37
ACCESSORIES	
ACCEDUNIED	41

PRELIMINARY EXAMINATION

Remove the instrument from the packing material and examine it carefully to make sure that no damage has occurred during shipping. If there is any damage, notify your Dealer or the nearest Hanna Customer Service Center. Each instrument is supplied with:

- HI 1131B Glass-body Combination pH Electrode with 1 m (3.3') Cable
- HI 7662 Temperature Probe
- HI 76404N Electrode Holder
- pH 4.01 & 7.01 Buffer Solutions (20 mL each)
- HI 7071 Electrolyte Solution
- HI 700661 Electrolyte Solution
- 12VDC Power Adapter
- Instruction Manual
- <u>Note</u>: Save all packing material until you are sure that the instrument functions correctly. All defective items must be returned in the original packing with the supplied accessories.

GENERAL DESCRIPTION

The HANNA **HI 2214** and **HI 2215** are a microprocessor-based pH and Temperature bench meter.

HI 2215 can also measure Oxidation Reduction Potential (ORP) in mV range. Relative mV feature is also provided.

pH measurements are compensated for temperature effect manually or automatically with the $HI\ 7662$ temperature probe.

Up to a five-point pH calibration can be performed using seven standard buffers. In addition, two custom buffers can be used during calibration.

The GLP feature provides data consistency.


A calibration due alarm can be set to alert the user that too much time elapsed since the last pH calibration.

Data can be stored in the meter's memory for later retrieval.


The meter's memory can hold a 100 manual logged points and 500 lot logging points (**HI 2215** only).

A USB connection ensures communication with a PC.

FUNCTIONAL DESCRIPTION



<u>Rear Panel</u>



9

1) Liquid Crystal Display (LCD).

(1) (2) (3) (4) (5) (6) (7) (8)

- 2) CAL key, to enter and exit calibration mode.
- 3) **CFM/GLP** key, to confirm calibration, different values or to display Good Laboratory Practice information.
- A °C key, to manually increase temperature value or other parameters.
- S) ♥ °C key, to manually decrease temperature value or other parameters.
- 6) **SETUP** key, to enter/exit SETUP mode.
- RANGE key, to select measurement range, switch to focused data in SETUP or toggle between buffer value and temperature during calibration.
- 8) LOG/CLR key, to store a value into memory, to clear pH calibration, or to delete log records
- 9) RCL key, memory recall.
- 10) **ON/OFF** switch.
- 11) Power supply socket.
- 12) USB connector.
- 13) BNC electrode connector.
- 14) Temperature probe socket.
- 15) Electrode reference socket.

HI 2214 SPECIFICATIONS

	—2.0 to 16.0 pH —2.00 to 16.00 pH
RANGE	±699.9 mV ±2000 mV
	-20.0 to 120.0°C (-4.0 to 248.0°F)
	0.1 pH, 0.01 pH
RESOLUTION	0.1 mV (±699.9mV) 1 mV (±2000mV)
	0.1 °C
	±0.01 pH
ACCURACY @ 20°C / 68°F	\pm 0.2 mV (\pm 699.9 mV) \pm 1 mV (\pm 2000 mV)
	± 0.2 °C (excluding probe error)
Rel mV offset range	± 2000 mV
pH Calibration	Automatic, up to 5 point calibration, 7 standard buffers available (1.68, 4.01, 6.86, 7.01, 9.18, 10.01, 12.45), and 2 custom buffers
Temperature compensation	Manual or Automatic from: —20.0 to 120.0 °C (—4.0 to 248.0 °F)
pH Electrode	HI 1131B (included)
Temperature probe	HI 7662 (included)
Input impedance	10 ¹² ohms
Log feature	100 records log on demand
PC communication	Opto-isolated USB
Power supply	12 VDC adapter (included)
Environment	0 — 50 °C (32 — 122 °F) max. 95% RH non-condensing
Dimensions	235 x 222 x 109 mm (9.2 x 8.7 x 4.3")
Weight	1.3 Kg (2.9 lbs); kit with holder 2.1 Kg (4.6 lb)
Warranty	2 years

HI 2215 SPECIFICATIONS

	—2.0 to 16.0 pH —2.00 to 16.00 pH —2.000 to 16.000 pH		
RANGE	\pm 999.9 mV (ISE & ORP) \pm 2000 mV (ISE & ORP)		
	—20.0 to 120.0 °C (-4.0 to 248.0°C)		
	0.1 pH, 0.01 pH, 0.001 pH		
RESOLUTION	0.1 mV (ISE & ORP) 1 mV (ISE & ORP)		
	0.1 °C		
	±0.01 pH		
ACCURACY @ 20°C / 68°F	±0.2 mV (±999.9 mV) ±1 mV (±2000 mV)		
	± 0.2 °C (excluding probe error)		
Rel mV offset range	$\pm 2000 \text{ mV}$		
pH Calibration	Automatic, up to 5 point calibration, 7 standard buffers available (1.68, 4.01, 6.86, 7.01, 9.18, 10.01, 12.45), and 2 custom buffers		
Temperature compensation	Manual or Automatic from: —20.0 to 120.0 °C (—4.0 to 248.0 °F)		
pH Electrode	HI 1131B glass-body, single junction, refillable, BNC conector (included)		
Temperature probe	HI 7662 (included)		
Input impedance	1012 ohms		
Log feature	100 records log on demand		
Log Interval feature	Stability logging (stab), 5, 10, 30 sec 1, 2, 5, 10, 15, 30, 60, 120, 180 min		
PC communication	Opto-isolated USB		
Power supply	12 VDC adapter (included)		
Environment	0 — 50 °C (32 — 122 °F) max. 95% RH non-condensing		
Dimensions	235 x 222 x 109 mm (9.2 x 8.7 x 4.3″)		
Weight	1.3 Kg (2.9 lbs); kit with holder 2.1 Kg (4.6 lb)		
Warranty	2 years		

OPERATIONAL GUIDE

POWER CONNECTION

Plug the 12 VDC adapter into the power supply socket.

- <u>Notes</u>: This instrument uses non volatile memory to retain the pH, mV, temperature calibrations and all other settings, even when unplugged.
 - Make sure a fuse protects the main line.

ELECTRODE AND PROBE CONNECTIONS

For pH or ORP measurements connect an electrode with internal reference to the BNC connector on the back of the instrument.

For electrodes with a separate reference connect the electrode's BNC to the BNC connector and the reference electrode plug to the reference socket. For temperature measurements and automatic temperature compensation connect the temperature probe to the appropriate socket.

INSTRUMENT START-UP

- Turn the instrument on by pressing the ON/OFF switch located on the rear panel.
- All LCD tags are displayed and a beep is sounded while the instruments perform a self test.



 The instrument will display "LoAd" message and "Z" blinking until initialization is complete.

pH MEASUREMENTS

Make sure the instrument has been calibrated before taking pH measurements.

- Submerse the electrode tip and the temperature probe approximately 3 cm (1¼") into the sample to be tested and stir gently. Allow time for the electrode to stabilize.
- The pH is displayed on the primary LCD and the temperature on the secondary LCD.



 If the reading is out of range, the closest full-scale value will be displayed blinking on the primary LCD. If measurements are taken successively in different samples, it is recommended to rinse the electrode thoroughly with deionized water or tap water and then with some of the next sample to prevent cross-contamination.

The pH reading is affected by temperature. In order to measure the pH accurately, the temperature effect must be compensated for. To use the Automatic Temperature Compensation feature, connect and submerse the **HI 7662** temperature probe into the sample as close as possible to the electrode and wait for a few seconds.

If the temperature of the sample is known, manual temperature compensation can be used by disconnecting the temperature probe.

The display will then show the last temperature reading with the "°C" tag blinking.

The temperature can now be adjusted with the **ARROW** keys (from -20.0 °C to 120.0 °C).





mV/ORP MEASUREMENTS

An optional ORP electrode must be used to perform ORP measurements (see Accessories).

Oxidation-Reduction Potential (REDOX) measurements provide the quantification of the oxidizing or reducing power of the tested sample.

The surface of the ORP electrode must be clean and smooth in order to obtain an accurate measurement.



(11/4"

- Press **RANGE** to enter mV range.
- Submerse the tip of the ORP electrode 3 cm (1¼") into the sample to be tested and allow a few seconds for the reading to stabilize.
- The instrument displays the mV reading on the primary LCD and the temperature on the secondary LCD.




• If the reading is out of range, the closest full-scale value will be displayed blinking on the primary LCD.



RELATIVE mV MEASUREMENTS

 Press RANGE until "rEL" message and "mV" tag are displayed for one second. After one second the temperature will be displayed on the secondary LCD and the "mV" tag will blink.



The reading displayed by the instrument is equal to the difference between the current mV input value and relative mV offset established in the relative mV calibration.

TEMPERATURE MEASUREMENTS

Connect the **HI 7662** temperature probe to the appropriate socket and turn the instrument on. Submerse the temperature probe into the sample and allow the reading on the secondary LCD to stabilize.



pH CALIBRATION

Calibrate the instrument frequently, especially if high accuracy is required.

The instrument should be recalibrated:

- Whenever the pH electrode is replaced.
- At least once a week.
- After testing aggressive chemicals.

• If "CAL" "INTV" tags are blinking during measurement.

Every time you calibrate the instrument use fresh buffers and perform an electrode Cleaning Procedure (see page 38).

PREPARATION

Pour small quantities of the buffer solutions into clean beakers. If possible, use plastic or glass beakers to minimize any EMC interferences. For accurate calibration and to minimize cross-contamination, use two beakers for each buffer solution. One for rinsing the electrode and one for calibration.

If you are measuring in the acidic range, use pH 7.01 or 6.86 as first buffer and pH 4.01 as second buffer. If you are measuring in the alkaline range, use pH 7.01 or 6.86 as first buffer and pH 10.01 or 9.18 as second buffer.

PROCEDURE

Calibration can be performed up to five-points.

For accurate measurements a three-point calibration is recommended. The calibration buffer can be selected from the calibration buffer list that includes custom buffers and the memorized standard buffers:

• pH 1.68, 4.01, 6.86, 7.01, 9.18, 10.01 and 12.45.

The custom buffers allow the user to calibrate in a buffer solution different from a standard one. Up to two custom buffers can be set in SETUP menu (see page 26). Each custom buffer value can be changed in a \pm 1.0 pH window around the set value, during calibration, when it is selected; the **"BUFFER pH"** tag will blink.

The instruments will automatically skip the buffer used during calibration and the buffers which are in a \pm 0.2 pH window, around one of the calibrated buffers.

All new calibrations will override existing stored calibration data in a \pm 0.2 pH window. The slopes adjacent to the new points will be reevaluated.

FIVE-POINT CALIBRATION

- Submerse the pH electrode and the temperature probe approximately 3 cm (1¼") into a buffer solution and stir gently. The temperature probe should be close to the pH electrode.
- Press CAL. The "CAL" and "
 "CAL" and "
 "1" tags will appear and the "7.01" buffer will be displayed on the secondary LCD.





- If necessary, press the ARROW keys to select a different buffer value.
 - The "Z" tag will blink on the LCD until the reading is stable.
- When the reading is stable and close to the selected buffer, the "READY" tag will be displayed and the "CFM" tag will blink.
- Press CFM to confirm calibration.

•



• The calibrated value will be displayed on the primary LCD and the second expected buffer value on the secondary LCD.



- After the first calibration point is confirmed, submerse the pH electrode and the temperature probe approximately 3 cm (1¼") into the second buffer solution and stir gently. The temperature probe should be close to the pH electrode.
- If necessary, press the **ARROW** keys to select a different buffer value.
- The "\Imes" tag will blink on the LCD until the reading is stable.
- When the reading is stable and close to the selected buffer, the "**READY**" tag will be displayed and the "**CFM**" tag will blink.



- Press CFM to confirm calibration.
- The calibrated value is then displayed on the primary LCD and the third expected buffer value on the secondary LCD.
- After the second calibration point is confirmed, submerse the pH electrode and the temperature probe approximately 3 cm (1¼") into the next buffer solution and stir gently. The temperature probe should be close to the pH electrode.
- If necessary, press the ARROW keys to select a different buffer value.
- The "\alpha" tag will blink on the LCD until the reading is stable.
- When the reading is stable and close to the selected buffer, the "READY" tag will be displayed and the "CFM" tag will blink.



CFM € READY

- Press **CFM** to confirm calibration.
- The calibrated value is then displayed on the primary LCD and the fourth expected value on the secondary LCD.
- After the third calibration point is confirmed, submerse the pH electrode and the temperature probe approximately 3 cm (1¼") into the next buffer solution and stir gently. The temperature probe should be close to the pH electrode.
- If necessary, press the ARROW keys to select a different buffer value.

- The "\alpha" tag will blink on the LCD until the reading is stable.
- When the reading is stable and close to the selected buffer, the "READY" tag will be displayed and the "CFM" tag will blink.
- Press CFM to confirm calibration.
- The calibrated value is then displayed on the primary LCD and the fifth expected value on the secondary LCD.
- After the fourth calibration point is confirmed, submerse the pH electrode and the temperature probe approximately 3 cm (11/4") into the next buffer solution and stir gently. The temperature probe should be close to the pH electrode.
- If necessary, press the ARROW keys to select a different buffer value.
- The "\alpha" tag will blink on the LCD until the reading is stable.
- When the reading is stable and close to the selected buffer, the "READY" tag will be displayed and the "CFM" tag will blink.
- Press CFM to confirm the fifth calibration point.
- The instrument stores the calibration value of the fifth point and returns to normal measurement mode.

FOUR, THREE OR TWO-POINT CALIBRATION

- Proceed as described in "FIVE-POINT CALIBRATION" section.
- Press CAL after the fourth, third or second calibration point was confirmed. The instrument will memorize the calibration data and return to measurement mode.

ONE-POINT CALIBRATION

Two SETUP selectable options are available: "Pnt" and "OFFS".

If the "**Pnt**" option is selected, the new calibration point overrides an existing one. The adjacent slopes will be reevaluated.

If the "**OFFS**" option is selected, an electrode offset correction is performed. The adjacent slopes will remain unchanged.

- Proceed as described in "FOUR, THREE or TWO-POINT CALIBRATION" section.
- Press CAL after the first calibration point was confirmed. The instrument will memorize the one-point calibration data and return to measurement mode.





- Notes: If the value measured by the instrument is not close to the selected buffer, "WRONG" "in and "WRONG" "I" tags will blink alternately. Check if the correct buffer has been used, or clean the electrode by following the Cleaning Procedure (see page 38). If necessary, change the buffer or the electrode.
 - When a custom buffer is displayed, the "BUFFER pH" tag blinks. To change the custom buffer value in accordance with the buffer temperature proceed as described in "WORKING WITH CUSTOM BUFFERS" (see page 14).
 - If the buffer temperature or the manual temperature exceeds the temperature limits of the buffer, "WRONG" tag and temperature reading will blink.
 - If "WRONG", "BUFFER pH" tags and "OLd" message are displayed blinking on the secondary LCD line, an inconsistency between new and previous (old) calibration is detected. Clear calibration parameters and proceed with calibration from the current calibration point. The instrument will keep all confirmed values during current calibration.
 - To clear calibration parameters for all buffers starting with current buffer, press CLR. The calibration will continue from the current point. If this procedure is performed while calibrating in the first point, the instrument returns to measurement mode.
 - Press RANGE to toggle between pH buffer, calibration buffer number and temperature reading.



• Each time a buffer is confirmed, the new calibration data replaces the old calibration data for the corresponding buffer. If current buffer has no previous data in the existing calibration and the calibration is not full (five buffers), the current buffer is added to the existing data. If the existing calibration is full, the instrument asks which buffer to be replace.

Press the **ARROW** keys to select another buffer to be replaced.

Press CFM to confirm the buffer that will be replaced.

Press CAL to leave calibration without replacing.

<u>Note</u>: If the replaced buffer is outside the \pm 0.2 pH window, around each of the calibrated buffers, it is possible to select this buffer for next calibration during current calibration.

WORKING WITH CUSTOM BUFFERS

If a custom buffer was set in SETUP menu, it can be selected during calibration by pressing the $\ensuremath{\mathsf{ARROW}}$ keys. The " $\ensuremath{\mathsf{BUFFER}}$ $\ensuremath{\mathsf{pH}}$ " tag will blink.

Press **SETUP** if you want to adjust the buffer value. The buffer value will start blinking.

Use the **ARROW** keys to change the buffer value.



After 5 seconds the buffer value is updated. Press **SETUP** if you want to change it again.

Note: Custom buffer value can be adjusted in a $\pm 1.00~\text{pH}$ window, around the set value.

pH BUFFER TEMPERATURE DEPENDENCE

The temperature has an effect on pH. The calibration buffer solutions are affected by temperature changes to a lesser degree than normal solutions. During calibration the instrument will automatically calibrate to the pH value corresponding to the measured or set temperature.

TE	MP	pH BUFFERS						
°C	٩F	1.68	4.01	6.86	7.01	9.18	10.01	12.45
0	32	1.67	4.01	6.98	7.13	9.46	10.32	13.38
5	41	1.67	4.00	6.95	7.10	9.39	10.24	13.18
10	50	1.67	4.00	6.92	7.07	9.33	10.18	12.99
15	59	1.67	4.00	6.90	7.05	9.27	10.12	12.80
20	68	1.68	4.00	6.88	7.03	9.22	10.06	12.62
25	77	1.68	4.01	6.86	7.01	9.18	10.01	12.45
30	86	1.68	4.02	6.85	7.00	9.14	9.96	12.29
35	95	1.69	4.03	6.84	6.99	9.11	9.92	12.13
40	104	1.69	4.04	6.84	6.98	9.07	9.88	11.98
45	113	1.70	4.05	6.83	6.98	9.04	9.85	11.83
50	122	1.71	4.06	6.83	6.98	9.01	9.82	11.70
55	131	1.72	4.08	6.84	6.98	8.99	9.79	11.57
60	140	1.72	4.09	6.84	6.98	8.97	9.77	11.44
65	149	1.73	4.11	6.84	6.99	8.95	9.76	11.32
70	158	1.74	4.12	6.85	6.99	8.93	9.75	11.21
75	167	1.76	4.14	6.86	7.00	8.91	9.74	11.10
80	176	1.77	4.16	6.87	7.01	8.89	9.74	11.00
85	185	1.78	4.17	6.87	7.02	8.87	9.74	10.91
90	194	1.79	4.19	6.88	7.03	8.85	9.75	10.82
95	203	1.81	4.20	6.89	7.04	8.83	9.76	10.73

During calibration the instrument will display the pH buffer value at 25 °C.

RELATIVE mV CALIBRATION

- Press CAL when the instrument is in RELATIVE mV measurement mode. The "mV" and "Z" tags will be displayed. Absolute mV is displayed on the primary LCD and "AbS" message is displayed on the secondary LCD.
- When the absolute reading is stable and in measurement range, the instrument asks for confirmation.
- If the reading is out of range, "WRONG" tag will be displayed.
- Press CFM to confirm the absolute value. The instrument will display "0.0 mV" on the primary LCD and "rEL" message on the secondary LCD. In this moment the relative mV offset is equal to absolute mV reading.
- Use the **ARROW** keys if you want to change the displayed relative mV value.
- Press CFM to confirm the relative mV value. The relative mV offset is displayed on the primary LCD. The instrument returns to measurement mode.
- <u>Note</u>: The relative mV value can be changed only inside the relative mV offset window (\pm 2000 mV).

GOOD LABORATORY PRACTICE (GLP)

GLP is a set of functions that allows storage and retrieval of data regarding the maintenance and status of the electrode.

All data regarding pH or Rel mV calibration is stored for the user to review when necessary.

EXPIRED CALIBRATION

These instruments allow the user to set the number of days before the next required calibration. This value can be set from 1 to 7 days. The default setting is off (disabled). The "CAL" "INTV" tags will start blinking to warn the user that the instrument should be recalibrated.

For example, if a 4 days time out has been selected, the instrument will issue the alarm exactly 4 days after the last calibration.

If the expiration value is changed (e.g. to 5 days), then the alarm will be immediately recalculated and appear 5 days after the last calibration.

- <u>Notes</u>: When the instrument is not calibrated or calibration is cleared (default values loaded) there is no "expired calibration", and the display always shows the "CAL" "INTV" tags blinking.
 - When an abnormal condition in the RTC is detected, the instrument forces the "expired calibration" status.

pH CALIBRATION DATA

Calibration data is stored automatically after a successful calibration.

To view the pH calibration data, press **GLP** while in measurement mode. The instruments will display the time (hh:mm) of the last calibration.




Use the **ARROW** keys to scroll through the calibration data:

• The date (yyyy.mm.dd).



• The pH calibration offset.



 The pH calibration slope (the GLP slope is the average of the calibration slopes; the percentage is referred to the ideal value of 59.16 mV/pH).



The pH calibration buffers in calibrating order and with the selected resolutions at calibration moment.

The first all a clibration to final terms to final terms to the first set of the selected terms to the first set of the selected terms to the selected terms terms to the selected terms terms terms to the selected terms term

The first pH calibration buffer:



The second pH calibration buffer:

The third pH calibration buffer:

The fourth pH calibration buffer:

The fifth pH calibration buffer:

- <u>Notes</u>: The "OLd" message displayed beside the pH value means that this buffer was not used during last calibration. Press and hold down SETUP if you want to see calibration date (or time, if old calibration was performed on the same day as the current calibration).
 - For each custom buffer used in calibration, the "BUFFER pH" tag will blink.
 - If "no bUF" message appears on the LCD, the instrument informs you that the calibration was performed with less than five buffers.



• Calibration Expiration status:

- if disabled.



- or the number of days until the calibration alarm will be displayed.



- or if expired (7 days ago).

• The instrument ID.

- Notes: Press GLP to return to measurement mode.
 - If calibration has not been performed, the instrument displays "no CAL" message blinking.



RELATIVE mV CALIBRATION DATA

Last Relative mV calibration data is stored automatically after a successful calibration.

To view the Relative mV calibration data, press **GLP** while in Relative mV measurement mode.

The instrument will display the Relative mV GLP information and the time (hh:mm:ss) of the last Rel mV calibration.

Use the **ARROW** keys to scroll through the calibration data:

- The Relative mV calibration date.
- The Relative mV calibration offset.
- The instrument ID.
- Notes: Press GLP to return to measurement mode.
 - If calibration has not been performed, the instrument displays "no CAL" message blinking.

LOGGING FUNCTION

Up to 600 logged samples can be stored into memory.

A 100 manually logged records and 500 lot logging records can be stored in the memory. To select logging type enter SETUP menu.


LOGGING THE CURRENT DATA (manual logging)

Select the manual logging mode in SETUP menu.

To store the current reading into memory press **LOG** while the instrument is in measurement mode.

The instrument will display:

HI 2214 - the current date (mm.dd) on the primary LCD, the record number on the secondary LCD and "**LOG**" tag will blink for a few seconds (see example below: record No. 15 dated September 14):




HI 2215 - "MAn" on the primary LCD, the record number on the secondary LCD and "LOG" tag for a few seconds (see example below: record No. 15):



followed by the number of free records:



If there are less than 6 memory locations remaining, the record number and "Lo" message will be displayed to alert the user.



If the log space is full (100 records), "FULL LOG" message will be displayed and no more data will be saved.



When **LOG** is pressed, a complete set of information is stored: date, time, pH, mV, temperature and calibration data.

LOT LOGGING (HI 2215 only)

Select "StAb" (for stability logging) or the desired time interval.

To start interval logging press ${\rm LOG}$ key while the instrument is in measurement mode.

When the selected interval is reached, the instrument will display the current lot number on the primary LCD line, the record number on the secondary LCD line and the "LOG" tag (see example below: Lot 5 record 7)



followed by the number of free records on the corresponding memory space.



If stability logging is selected, a complete set of data is memorized every time the reading becomes stable after an unstable condition.

To stop interval logging press LOG key again.

<u>Note</u>: When pressing any key that is not active, while Log interval is running, the following message is displayed for a few seconds:



VIEW LOGGED DATA

Press **RCL** while in measurement mode to retrive the stored information. If no data was logged for the selected measurement range and no lots are memorized (**HI 2215** only), one of the next messages will be displayed:

No records:



No pH measurement records:



No Relative mV and mV records:



Otherwise:

- For HI 2214 the instrument will display the pH or relative mV value on the primary LCD and the last stored record number on the secondary LCD, along with "LOG" tag.
- For HI 2215 the instrument will display the lot number on the primary LCD, the number of records on the secondary LCD, "LOG" tag and "CFM" blinking. If samples were logged on demand "MAn" will be displayed on the primary LCD and the number of samples logged on the secondary LCD (see example below: manual log, 15 samples logged).



Press **ARROW** keys to select different lot.



All the existing lots are displayed.

Press CFM to view record information.

The instrument will display the **pH** or **relative mV** value on the primary LCD and the last stored record number on the secondary LCD, along with "LOG" tag.




Oľ

LOG 352.7	ЧЧ
-----------	----

Use the **ARROW** keys to scroll through the records. To view additional information press **RANGE**:

• The **mV** value is displayed on the primary LCD and the temperature value on the secondary LCD.




• The time (hh.mm.ss) is displayed on the LCD, along with the "TIME" tag.



 The date (yyyy.mm.dd) is displayed on the LCD along with the "DATE" tag.



 The calibration offset is displayed on the primary LCD and "OFS" message on the secondary LCD.



 The calibration slope is displayed on the primary LCD and "SLP" message on the secondary LCD.



- <u>Note</u>: If in Relative mV RECALL mode, the instruments will display "----" message on the primary LCD.
- The interval for lot logging (HI 2215 only).




To delete manual logged records press **CFM** while "**MAn**" is displayed to view manually logged records. Press the **CLR** key, "**dEL**" and the record number will be displayed. Press **CFM** to delete. Use the **ARROW** keys to change the record number.

To delete a lot use the **ARROW** keys to select the desired lot. Press **CLR** key, "**dEL Lot**" will appear on the display. Press **CFM** key to delete.

The "dEL" message is displayed on the primary LCD and the selected record on the secondary LCD, along with "LOG" tag.



- Notes: The ARROW keys can be used to change the record.
 - Press SETUP to delete all records/lots. The display will show "dEL" in the primary LCD and "ALL" in the secondary LCD.



- Press CFM to confirm delete.
- Press CAL or RANGE to escape and return to RCL screen.
- If "dEL ALL" option was selected, all the log on demand records or lots are deleted. While deleting the "Σ" tag is displayed blinking.
- Press RCL to exit record information and enter lot information (HI 2215 only).
- Press RCL again to return to measurement mode.
- If one or more records/lots were deleted the "\"a" tag blinks until the log memory space is reorganized.

SETUP

Setup mode allows viewing and modifying the following parameters:

- Expired Calibration Alarm
- Log Interval (HI 2215 only)
- First Custom Buffer
- Second Custom Buffer
- One-point Calibration Behavior
- pH Resolution
- Current Time (hour & minute)
- Current Date (year, month & day)
- Beep Status
- Instrument Id
- Temperature Unit

To enter SETUP mode press **SETUP** while the instrument is in measurement mode.

Select a parameter with the **ARROW** keys.

Press **CAL** to change a parameter value. The selected parameter will start blinking.

Press **RANGE** to toggle between displayed parameters.

Press the **ARROW** keys to increase or decrease the displayed value.

Press CFM to save the modified value or CAL to escape without saving.


EXPIRED CALIBRATION ALARM

Press CAL when the calibration time-out is displayed. Calibration time-out ("OFF" or "1" to "7" days) and "CFM" will start blinking.

Press the **ARROW** keys to change the calibration time-out value.

Press **CFM** to save the modified calibration time-out value.

Press CAL to escape without saving.



LOG INTERVAL (HI 2215 only)

Press **CAL** when log interval is displayed. The log interval and "**CFM**" will start blinking ("**MAn**" for log on demand, "**StAb**" log on stability, interval in seconds or minutes).



Press the **ARROW** keys to change the custom buffer value. Press **CFM** to confirm the selection. Press **CAL** to escape without saving.

FIRST CUSTOM BUFFER

Press **CAL** when "**cb1**" is displayed. The custom buffer (disabled — "**no**" or "**0**" to "**16**" pH) and "**CFM**" will start blinking.

CAL	pН	
ECFM =		
		i.

Press the **ARROW** keys to change the custom buffer value.


Press CFM to save the modified custom buffer value.

Press CAL to escape without saving.

 $\underline{\textit{Note}}$: To delete custom buffer enter change mode by pressing CAL then

press CLR key. SECOND CUSTOM BUFFER

Press **CAL** when "**cb2**" is displayed. The custom buffer (disabled – "**no**" or "**0**" to "**16**" pH) and "**CFM**" will start blinking.



Press the **ARROW** keys to change the custom buffer value.

Press CFM to save the modified custom buffer value.

Press **CAL** to escape without saving.

ONE-POINT CALIBRATION BEHAVIOR

Press CAL when "1 Pnt" message is displayed on the secondary LCD. One of the two options ("Pnt" or "OFFS") and "CFM" will start blinking (see pH CALIBRATION PROCEDURE for details, page 9).



Press the ARROW keys to toggle between "Pnt" and "OFFS" options.

Press **CFM** to save the behavior for one-point calibration.

Press **CAL** to escape without saving.

pH RESOLUTION

Press **CAL** when "**rES**" message is displayed on the secondary LCD. The set resolution ("0.1", "0.01" or "0.001") and "**CFM**" will start blinking.


Press the **ARROW** keys to toggle between 0.1, 0.01 and 0.001 options for **HI 2215**; 0.1 and 0.01 options for **HI 2214**.

Press **CFM** to save the modified value.

Press CAL to escape without saving.


CURRENT TIME

Press **CAL** when the current time is displayed. The hour and "**CFM**" will start blinking.



Press the **ARROW** keys to change the hour.

Press RANGE. The minutes and "CFM" will start blinking.



Press the **ARROW** keys to change the minutes.

Press **CFM** to save the modified value.

Press **CAL** to escape without saving.

CURRENT DATE

Press ${\bf CAL}$ when the current date is displayed. The year and ${\bf ``CFM''}$ will start blinking.



Press the **ARROW** keys to change the year.

Press RANGE. The month and "CFM" will start blinking.



Press the **ARROW** keys to change the month. Press **RANGE**. The day and "**CFM**" will start blinking.



Press the **ARROW** keys to change the day. Press **CFM** to save the modified value.

Press CAL to escape without saving.

Tess CAL to escupe without saving.

BEEP STATUS

Press **CAL** when the beep status is displayed. Beep status ("**On**" or "**OFF**") and "**CFM**" will start blinking.



Press the ARROW keys to change the beep status (On or OFF).

Press **CFM** to save the modified beep status.

Press CAL to escape without saving.

When enable, beep sounds as a short beep every time a key is pressed or when the calibration can be confirmed.

A long beep alert that the pressed key is not active or a wrong condition is detected while in calibration.

INSTRUMENT ID

Press **CAL** when the **"InId**" is displayed. The instrument ID (**"0000"** to **"9999"**) and **"CFM**" will start blinking.



Press the **ARROW** keys to change the instrument ID value.

Press CFM to save the modified instrument ID value.

Press CAL to escape without saving.

<u>Note</u>: The instrument ID is downloaded to a PC as part of a logged data, set to identify its origin.

TEMPERATURE UNIT

Press CAL when "tnP" is displayed. The temperature unit and "CFM" will start blinking.



Press the **ARROW** keys to change the option. Press **CFM** to save the modified temperature unit. Press **CAL** to escape without saving.

TEMPERATURE CALIBRATION (for technical personnel only)

All the instruments are factory calibrated for temperature.

Hanna's temperature probes are interchangeable and no temperature calibration is needed when they are replaced.

If the temperature measurements are inaccurate, temperature recalibration should be performed.

For an accurate recalibration, contact your dealer or the nearest Hanna Customer Service Center, or follow the instructions below.

- Prepare a vessel containing ice and water and another one containing hot water (around 50 °C). Place insulation material around the vessels to minimize temperature changes.
- Use a calibrated thermometer with a resolution of 0.1°C as a reference thermometer.
- With the instrument off, press and hold down the CFM & SETUP keys, then power on the instrument. The "CAL" tag will appear and the secondary LCD will show "0.0 °C".



- Submerse the temperature probe in the vessel with ice and water as near as possible to the reference thermometer. Allow a few seconds for the probe to stabilize.
- Use the ARROW keys to set the reading on the secondary LCD to that
 of ice and water, measured by the reference thermometer. When the
 reading is stable and close to the selected calibration point, "READY"
 tag will appear and "CFM" tag will blink.
- Press CFM to confirm. The secondary LCD will show "50.0 °C".



 Submerse the temperature probe in the second vessel as near as possible to the reference thermometer. Allow a few seconds for the probe to stabilize. • Use the **ARROW** keys to set the reading on the secondary LCD to that of the hot water.



- When the reading is stable and close to the selected calibration point, "READY" tag will appear and "CFM" tag will blink.
- Press CFM to confirm. The instrument memorize calibration and restart to measurement mode.
- <u>Note</u>: If the reading is not close to the selected calibration point, "WRONG" tag will blink. Change the temperature probe and restart calibration.

mV CALIBRATION (for technical personnel only)

All the instruments are factory calibrated for mV.

Hanna's ORP electrodes are interchangeable and no mV calibration is needed when they are replaced.

If the mV measurements are inaccurate, mV recalibration should be performed.

For an accurate recalibration, contact your dealer or the nearest Hanna Customer Service Center, or follow the instructions below.

A two-point calibration can be performed at 0.0 mV and 1800.0 mV.

- Attach to the BNC connector a mV simulator with an accuracy of ± 0.1 mV.
- With the instrument off, press and hold down the CAL & ▼ °C keys, then power on the instrument. The "CAL" tag will appear and the secondary LCD will show "0.0 mV".
- Set 0.0 mV on the simulator. When the reading is stable and close to the selected calibration point, "READY" tag will appear and "CFM" tag will blink.
- Press CFM to confirm. The secondary LCD will display "1800 mV".
- Set 1800.0 mV on the simulator. When the reading is stable and close to the selected calibration point, "READY" tag will appear and "CFM" tag will blink.
- Press CFM to confirm. The instrument memorize calibration and restart to measurement mode.

- <u>Notes</u>: If the reading is not close to the selected calibration point, "WRONG" tag will blink. Verify calibration condition or contact your vendor if you can not calibrate.
 - Pressing CAL key during temperature or mV calibration the instrument quit calibration mode and restart to measurement mode without memorizing calibration.

PC INTERFACE

Data transmission from the instrument to the PC can be done with the HI 92000 Windows[®] compatible software (optional). HI 92000 also offers graphing and an on-line help feature.

Data can be exported to the most popular spreadsheet programs for further analysis.

To connect your instrument to a PC, use a standard USB cable connector. Make sure that your instrument is switched off and plug one connector to the instrument's USB socket and the other to the USB port of your PC.

<u>Note</u>: If you are not using Hanna Instruments **HI 92000** software, please see the following instructions.

SENDING COMMANDS FROM PC

It is also possible to remotely control the instrument with any terminal program. Use a standard USB cable to connect the instrument to a PC, start the terminal program and set the communication options as follows: 8, N, 1, no flow control, 9600 baud rate.

COMMAND TYPES

To send a command to the instrument follow the next scheme:

	<command prefix $>$ $<$ command $>$ $<$ CR $>$
where:	<command prefix $>$ is a 16 ASCII character
	<command $>$ is the command code.

Note: Either small or capital letters can be used.

SIMPLE COMMANDS

RNG	Is equivalent to pressing RANGE
CAL	Is equivalent to pressing CAL
CFM	Is equivalent to pressing CFM
UPC	Is equivalent to pressing the UP arrow key
DWC	Is equivalent to pressing the DOWN arrow key
LOG	Is equivalent to pressing LOG
RCL	Is equivalent to pressing RCL
SET	Is equivalent to pressing SETUP

- CHR xx Change the instrument range according with the parameter value (xx):
 - xx=00 pH range/0.001 resolution
 - xx=01 pH range/0.01 resolution
 - xx=02 pH range/0.1 resolution
 - xx=03 mV range
 - xx=04 Relative mV range

The instrument will answer for these commands with:

<STX> <answer> <ETX>

- where: <STX> is 02 ASCII code character (start of text)
 - <ETX> is 03 ASCII code character (end of text)

<answer>:

- <ACK> is 06 ASCII code character (recognized command)
- <NAK> is 21 ASCII code character (unrecognized command)
- <CAN> is 24 ASCII code character (corrupted command)

COMMANDS REQUIRING AN ANSWER

The instrument will answer for these commands with:

<STX> <answer> <checksum> <ETX>

where the checksum is the bytes sum of the answer string sent as 2 ASCII characters.

All the answer messages are with ASCII characters.

- **RAS** Causes the instrument to send a complete set of readings in according with the current range:
 - pH, temperature and mV reading on pH range.
 - Rel mV, absolute mV and temperature reading on Relative mV range.

The answer string contains:

- Meter mode (2 chars):
 - 00 pH range (0.001 resolution)
 - 01 pH range (0.01 resolution)
 - 02 pH range (0.1 resolution)
 - 03 mV ranae
 - 04 Rel mV range
- Meter status (2 chars of status byte): represents a 8 bit hexadecimal encoding.
 - 0x10 temperature probe is connected
 - 0x01 new GLP data available
 - 0x02 new SETUP parameter

- Reading status (2 chars): R in range, O over range, U - under range. First character corresponds to the appropriate range reading. Second character corresponds to mV reading.
- Primary reading (corresponding to the selected range)
 7 ASCII chars, including sign and decimal point.
- Secondary reading (only when primary reading is not mV) - 7 ASCII chars, including sign and decimal point.
- Temperature reading 8 ASCII chars, with sign and two decimal points, always in °C.
- MDR Requests the instrument model name and firmware code (16 ASCII chars).
- GLP Requests the calibration data record.

The answer string contains:

- GLP status (1 char): represents a 4 bit hexadecimal encoding.
 - 0x01 pH calibration available
 - 0x02 Rel mV calibration available
- pH calibration data (if available), which contains:
 - the number of calibrated buffers (1 char)
 - the offset, with sign and decimal point (7 chars)
 - the average of slopes, with sign and decimal point (7 chars)
 - the calibration time, yymmddhhmmss (12 chars)
 - buffers information (for each buffer)
 - type (1 char): 0 standard, 1 custom
 - status (1 char): N (new) calibrated in last calibration; O (old) - from an old calibration.
 - warnings during calibration (2 chars): 00 no warning, 04 - Clean Electrode warning.
 - buffer value, with sign and decimal point (7 chars).
 - the calibration time, yymmddhhmmss (12 chars).
 - electrode condition, with sign (3 chars). The "-01" code means not calculated.
- Rel mV calibration data (if available), which contains:
 - the calibration offset, with sign (7 chars)
 - the calibration time, yymmddhhmmss (12 chars).
- Requests the setup parameters setting.

The answer string contains:

PAR

- Instrument ID (4 chars)
- Calibration alarm time out for pH range (2 chars)
- SETUP information (2 chars): 8 bit hexadecimal encoding.
 - 0x01 beep ON (else OFF)
 - 0x04 degrees Celsius (else degrees Fahrenheit)
 - 0x08 Offset calibration (else Point calibration)
- The number of custom buffers (1 char)

Requests the number of logged samples (4 chars).

- The custom buffer values, with sign and decimal point, for each defined custom buffer (7 chars)
- NSLx
- Command Parameter (1 char):
 - x range (see Note)

LODPxxx Requests the xxxth pH record logged data.

LODMxxx Requests the xxxth mV/Rel mV record logged data.

LODPALLff Requests all pH Log on demand - ff=frame number

LODMALLff Requests all mV/Rel mV Log on demand - ff = frame number

- <u>Note</u>: The log on demand can be downloaded in frames of 8 records (01 is first frame number)
- LLSxff: Requests information about all lots in frames of 10 lots (a frame contains information about 10 lots) (HI 2215 only) Command Parameters:
 - x range (see Note)
 - ff requested frame number first frame is labeled 01
- GLDxxxff: Requests the records of the "xxx" lot number (HI 2215 only). The records are sent in frames of 10 records; "ff" is the frame number (01 first frame). (Example: Lot 13 has 53 records. The records will be sent in 6 frame, 5 with10 records and 1 with 3 records.)


Command Parameters:

- xxx Lot number (eq: for lot number 1 xxx = 001)
- ff requested frame number first frame is labeled 01

Errors: • "Err3" is sent if the Log on demand is empty.

- "Err4" is sent if the requested set parameter is not available.
- "Err5" is sent if an argument of the command is not correct.
- "Err6" is sent if the requested range is not available.
- "Err7" is sent if the instrument is in logging mode.
- "Err8" is sent if the instrument is not in measurement mode.
- Invalid commands will be ignored.
- Note: P request for pH range.
 - M request for mV and Rel mV ranges.

ELECTRODE CONDITIONING & MAINTENANCE



PREPARATION PROCEDURE

Remove the protective cap of the pH electrode.

DO NOT BE ALARMED IF SALT DEPOSITS ARE PRESENT. This is normal with electrodes. They will disappear when rinsed with water.

During transport, tiny bubbles of air may form inside the glass bulb affecting proper functioning of the electrode. These bubbles can be removed by "shaking down" the electrode as you would do with a glass thermometer.

If the bulb and/or junction is dry, soak the electrode in **HI 70300** or **HI 80300** Storage Solution for at least one hour.

For refillable electrodes:

If the filling solution (electrolyte) is more than 2½ cm (1") below the fill hole, add **HI 7082** or **HI 8082** 3.5M KCI Electrolyte Solution for double junction or **HI 7071** or **HI 8071** 3.5M KCI + AgCI Electrolyte Solution for single junction electrodes.

For faster response, unscrew the fill hole screw during measurements.

For AMPHEL® electrodes:

If the electrode does not respond to pH changes, the battery is run down and the electrode should be replaced.

MEASUREMENT

Rinse the electrode tip with distilled water. Submerse the tip (3 cm $/1 \frac{1}{4''}$) in the sample and stir gently for a few seconds.

For a faster response and to avoid cross-contamination of the samples, rinse the electrode tip with a few drops of the solution to be tested, before taking measurements.

STORAGE PROCEDURE

To minimize clogging and assure a quick response time, the glass bulb and the junction should be kept moist and not allowed to dry out.

Replace the solution in the protective cap with a few drops of **HI 70300** or **HI 80300** Storage Solution or, in its absence, Filling Solution (**HI 7071** or **HI 8071** for single junction and **HI 7082** or **HI 8082** for double junction electrodes). Follow the Preparation Procedure before taking measurements.

Note: NEVER STORE THE ELECTRODE IN DISTILLED OR DEIONIZED WATER.

PERIODIC MAINTENANCE

Inspect the electrode and the cable. The cable used for connection to the instrument must be intact and there must be no points of broken insulation on the cable or cracks on the electrode stem or bulb. Connectors must be perfectly clean and dry. If any scratches or cracks are present, replace the electrode. Rinse off any salt deposits with water.

For refillable electrodes:

Refill the reference chamber with fresh electrolyte (**HI 7071** or **HI 8071** for single junction and **HI 7082** or **HI 8082** for double junction electrodes). Allow the electrode to stand upright for 1 hour. Follow the Storage Procedure above.

CLEANING PROCEDURE

- General Soak in Hanna **HI 7061** or **HI 8061** General Cleaning Solution for approximately ½ hour.
- Protein Soak in Hanna HI 7073 or HI 8073 Protein Cleaning Solution for 15 minutes.
- Inorganic Soak in Hanna **HI 7074** Inorganic Cleaning Solution for 15 minutes.
- Oil/grease Rinse with Hanna **HI 7077** or **HI 8077** Oil and Fat Cleaning Solution.


IMPORTANT: After performing any of the cleaning procedures, rinse the electrode thoroughly with distilled water, refill the reference chamber with fresh electrolyte (not necessary for gel-filled electrodes) and soak the electrode in **HI 70300** or **HI 80300** Storage Solution for at least 1 hour before taking measurements.

TROUBLESHOOTING GUIDE

SYMPTOMS	PROBLEM	SOLUTION
Slow response/excessive drift.	Dirty pH electrode.	Clean the electrode and then soak the tip in HI 7061 or HI 8061 for 30 minutes.
Readings fluctuate up and down (noise).	Clogged/dirty junction. Low electrolyte level (refillable electrodes only).	Clean the electrode. Refill with fresh solution (for refillable electrodes only). Check cable and connectors.
The meter does not accept the buffer solution for calibration.	Dirty electrode or contaminated buffer.	Follow the cleaning procedure. If still no results, replace the electrode. Replace buffer.
If the display shows: "pH" and "-2.00" or "16.00" blinking.	Out of range in the pH scale.	 a) Verify that the electrode is connected. b) Verify that the shipping cap has been removed. c) Recalibrate the meter. d) Make sure the pH sample is in the specified range. e) Check electrolyte level and general state of the electrode.
If the display shows: "mV" and "-2000" or "2000" blinking	Out of range in the mV scale.	Verify that the electrode is connected.
The meter does not work with the temperature probe.	Broken temperature probe. Wrong temperature probe used.	Replace the temperature probe.
The meter fails to calibrate or gives faulty readings.	Broken pH electrode.	Replace the electrode.
At startup the meter displays all LCD tags permanently.	One of the keys is stuck.	Check the keyboard or contact the vendor.
" Err xx " error message displayed.	Internal error.	Power off the meter and then power it on. If the error persists, contact the vendor.

TEMPERATURE CORRELATION FOR pH SENSITIVE GLASS

The resistance of glass electrodes partially depends on the temperature. The lower the temperature, the higher the resistance. It takes more time for the reading to stabilize if the resistance is higher. In addition, the response time will suffer to a greater degree at temperatures below 25 °C.



Since the resistance of the pH electrode is in the range of 50 - 200 Mohm, the current across the membrane is in the pico Ampere range. Large currents can disturb the calibration of the electrode for many hours.

For these reasons high humidity environments, short circuits and static discharges are detrimental to a stable pH reading.

The pH electrode's life also depends on the temperature. If constantly used at high temperatures, the electrode life is drastically reduced. **Typical Electrode Life**

Ambient Temperature	1 — 3 years
90 °C	Less than 4 months
120 °C	Less than 1 month

Alkaline Error

High concentrations of sodium ions interfere with readings in alkaline solutions. The pH at which the interference starts to be significant depends upon the composition of the glass. This interference is called alkaline error and causes the pH to be underestimated. Hanna's glass formulations have the indicated characteristics.

ionnoidhean annoidhean annoidhean annoidhean annoidhean ann ann ann ann ann ann ann ann ann		
Sodium Ion Correction for the Glass at 20-25 °C		
Concentration	pH	Error
0.1 Mol L ⁻¹ Na+	13.00	0.10
	13.50	0.14
	14.00	0.20
	12.50	0.10
	13.00	0.18
1.0 Mol L-1 Na +	13.50	0.29
	14.00	0.40

ACCESSORIES

pH BUFFER SOLUTIONS

HI 70004P pH 4.01 Buffer Sachets, 20 mL, 25 pcs

- HI 70007P pH 7.01 Buffer Sachets, 20 mL, 25 pcs
- HI 70010P pH 10.01 Buffer Sachets, 20 mL, 25 pcs
- HI 7001L pH 1.68 Buffer Solution, 500 mL
- HI 7004L pH 4.01 Buffer Solution, 500 mL
- HI 7006L pH 6.86 Buffer Solution, 500 mL
- HI 7007L pH 7.01 Buffer Solution, 500 mL
- HI 7009L pH 9.18 Buffer Solution, 500 mL
- HI 7010L pH 10.01 Buffer Solution, 500 mL
- HI 8004L pH 4.01 Buffer Solution in FDA approved bottle, 500 mL
- HI 8006L pH 6.86 Buffer Solution in FDA approved bottle, 500 mL
- HI 8007L pH 7.01 Buffer Solution in FDA approved bottle, 500 mL
- HI 8009L pH 9.18 Buffer Solution in FDA approved bottle, 500 mL
- HI 8010L pH 10.01 Buffer Solution in FDA approved bottle, 500 mL

ELECTRODE STORAGE SOLUTIONS

- HI 70300L Storage Solution, 500 mL
- HI 80300L Storage Solution in FDA approved bottle, 500 mL

ELECTRODE CLEANING SOLUTIONS

- HI 70000P Electrode Rinse Sachets, 20 mL, 25 pcs
- HI 7061L General Cleaning Solution, 500 mL
- HI 7073L Protein Cleaning Solution, 500 mL
- HI 7074L Inorganic Cleaning Solution, 500 mL
- HI 7077L Oil & Fat Cleaning Solution, 500 mL
- HI 8061L General Cleaning Solution in FDA approved bottle, 500 mL
- HI 8073L Protein Cleaning Solution in FDA approved bottle, 500 mL
- HI 8077L Oil & Fat Cleaning Solution in FDA approved bottle, 500 mL

ELECTRODE REFILL ELECTROLYTE SOLUTIONS

- HI 7071 3.5M KCl + AgCl Electrolyte, 4x30 mL, for single junction electrodes
- HI 7072 1M KNO, Electrolyte, 4x30 mL
- HI 7082 3.5M KCl Electrolyte, 4x30 mL, for double junction electrodes
- HI 8071 3.5M KCl + AgCl Electrolyte in FDA approved bottle, 4x30 mL, for single junction electrodes
- HI 8072 1M KNO₂ Electrolyte in FDA approved bottle, 4x30 mL
- HI 8082 3.5M KCI Electrolyte in FDA approved bottle, 4x30 mL, for double junction electrodes.

ORP PRETREATMENT SOLUTIONS

HI 7091L	Reducing Pretreatmen	it Solution, 500 ml	
HI 7092L	Oxidizing Pretreatmen	it Solution, 500 ml	L


pH ELECTRODES

All electrodes part numbers ending in B are supplied with a BNC connector and 1 m (3.3') cable, as shown below :




HI 1043B

Glass-body, double junction, refillable, combination **pH** electrode. Use: strong acid/alkali.




HI 1053B

Glass-body, triple ceramic, conic shape, refillable, combination $\ensuremath{\textbf{pH}}$ electrode. Use: emulsions.




HI 1083B

Glass-body, micro, Viscolene, non-refillable, combination **pH** electrode. Use: biotechnology, micro titration.




HI 1131B

Glass-body, single junction, refillable, combination **pH** electrode. Use: general purpose.




HI 1330B

Glass-body, semimicro, single junction, refillable, combination **pH** electrode. Use: laboratory, vials.




HI 1331B

Glass-body, semimicro, single junction, refillable, combination **pH** electrode. Use: flasks.




HI 1230B

Plastic-body (**PES**), double junction, gel-filled, combination **pH** electrode. Use: general, field.



HI 2031B

Glass-body, semimicro, conic, refillable, combination **pH** electrode. Use: semisolid products.




HI 1332B

Plastic-body (**PES**), double junction, refillable, combination **pH** electrode. Use: general purpose.




FC 100B

Plastic-body (**PVDF**), double junction, refillable, combination **pH** electrode. Use: general purpose for food industry.




FC 200B

Plastic-body (PVDF), open junction, conic, Viscolene, non-refillable, combination **pH** electrode. Use: meat & cheese.




FC 210B

Glass-body, double junction, conic, Viscolene, non-refillable, combination ${\bf p}{\bf H}$ electrode. Use: milk, yagurt.




FC 220B

Glass-body, triple-ceramic, single junction, refillable, combination **pH** electrode. Use: food processing.




FC 911B

Plastic-body (**PVDF**), double junction, refillable with built-in amplifier, combination **pH** electrode. Use: very high humidity.



HI 1413B


Glass-body, single junction, flat tip, Viscolene, non-refillable, combination ${\bf p}{\bf H}$ electrode. Use: surface measurement.



ORP ELECTRODES


HI 3131B

Glass-body, refillable, combination platinum **ORP** electrode. Use: titration.




HI 3230B

Plastic-body (PES), gel-filled, combination platinum ORP electrode. Use: general purpose.



HI 4430B

Plastic-body (**PES**), gel-filled, combination gold **ORP** electrode. Use: general purpose.



Consult the Hanna General Catalog for more electrodes with screw-type or BNC connectors.

EXTENSION CABLE FOR SCREW-TYPE ELECTRODES (SCREW TO BNC ADAPTER)


HI 7855/1 Extension cable 1 m (3.3') long

HI 7855/3 Extension cable 3 m (9.9') long

	HI 7855 SERIES CABLE CONNECTORS CONNECTOR AND 3.0 mm (0.12") CABLE WITH BI	NC
CONNECT TO SCREW TYPE ELECTRODES		CONNECT TO THE BNC SOCKET OF THE METER

OTHER ACCESSORIES

- HI 710005 Voltage adapter from 115 VAC to 12 VDC (USA plug)
- HI 710006 Voltage adapter from 230 VAC to 12 VDC (European plug)
- HI 710012 Voltage adapter from 240 VAC to 12 VDC (UK plug)
- HI 710013 Voltage adapter from 230 VAC to 12 VDC (South Africa plug)
- HI 710014Voltage adapter from 230 VAC to 12 VDC (Australia plug)HI 76404NElectrode holder



HI 8427	pH and ORP electrode simulator with 1 m (3.3') coaxial cable ending in female BNC connectors
HI 931001	pH and ORP electrode simulator with LCD and 1 m (3.3') coaxial cable ending in female BNC connectors
HI 7662 HI 92000	Temperature probe with 1 m (3.3') cable Windows® compatible software.

Windows® is registered Trademark of "Microsoft Co."

RECOMMENDATIONS FOR USERS

Before using these products, make sure they are entirely suitable for the environment in which they are used.

Operation of these instruments in residential areas could cause unacceptable interferences to radio and TV equipment, requiring the operator to follow all necessary steps to correct interferences.

The glass bulb at the end of the pH electrode is sensitive to electrostatic discharges. Avoid touching this glass bulb at all times.

During operation, ESD wrist straps should be worn to avoid possible damage to the electrode by electrostatic discharges.

Any variation introduced by the user to the supplied equipment may degrade the instruments' EMC performance.

To avoid electrical shock, do not use these instruments when voltages at the measurement surface exceed 24 VAC or 60 VDC.

To avoid damage or burns, do not perform any measurement in microwave ovens.

Hanna Instruments reserves the right to modify the design, construction and appearance of its products without advance notice.



Hanna Instruments Inc.

Highland Industrial Park 584 Park East Drive Woonsocket, RI 02895 USA

Technical Support for Customers

Tel. (800) 426 6287 Fax (401) 765 7575 E-mail tech@hannainst.com www.hannainst.com

Local Sales and Customer Service Office

Printed in EUROPE (ROMANIA)

MAN2215 02/11