

Euroseal®

A simple, versatile thermosealing machine



Euroseal® thermosealing machine is attractive, compact, ergonomic and easy to use. Thanks to the transparent polycarbonate cover, sealing operations can be performed with even greater precision. The height and depth of the special, patented roll holder can be adjusted and it can also be wall-mounted, in order to meet the specific needs of each dental practice. The integrated support surface aids the sealing process and provides more practical working conditions. Anti-roll return blocking system. Safety system to prevent the roll from burning. Sound and light warning signals during use. Fully compliant with the UNI 868-5 and EN 11607-2 standards.



Sealing strip:

12 mm

Sealing area width:

310 mm



Dimensions (wxhxd):

Standard:
440 x 390 x 315 mm
Mounted roll holder:
440 x 300 x 200 mm
Double roll holder:
440 x 390 x 420 mm



Weight:

7.1 Kg



Power supply:

200 - 240 V
50/60 Hz

Absorbed power:

100 W



Versatile

The patented roll holder is adjustable and it can be wall-mounted.



Resistant

The light, cutting-edge materials can withstand extreme wear and tear.



Ergonomic

Every detail is designed to improve the quality of the operators' work.



Simple

Thanks to the transparent polycarbonate cover, rolls can be monitored down to the sealing area and below.



Practical

Working with the machine is easy, practical and handy thanks to the spacious support surface.



Safe

There are burning prevention and anti-roll return blocking systems, as well as sound and light signals.