



Labworldsoft® 6 Advanced

/// Технический паспорт

Labworldsoft® 6 Advanced lab software is unique on the market and the cockpit for your digital laboratory. On the PC, you are able to automate laboratory processes, store your own recipes and process flows or perform calculations during the experiment.

labworldsoft® 6 Advanced at a glance

- > Suitable for all device configurations
- > Access to more than 400 devices from the device library
- > Any number of devices can be controlled at the



designed for scientists

same time

- > For laboratory equipment from all manufacturers
- > 1 year of support with new versions of labworldsoft® Advanced

Your benefits

- > Time savings and increased safety through automation of processes
- > Great flexibility through individual configuration of laboratory processes and set-ups with devices from different manufacturers
- > Convenient user interface for a clear display of interlinked work processes
- > Easy user management

Try it free of charge

You can test labworldsoft® 6 Advanced for 30 days free of charge: During this time, you have access to the entire device library with over 400 laboratory devices as well as to all functions that are required for your new digital laboratory. There is no obligation for you. After 30 days, you either choose one of our licenses or your access will expire.



designed for scientists

Технические данные

Unlimited access to the device library	да
IKA Rotavisc	да
Circulators	да
pH-measurement devices	да
IKA Screening System	да
Rohde	

