



Labworldsoft® 6 Visc

/// Технический паспорт

Labworldsoft 6 Visc software enables the automation of laboratory processes by connecting ROTAVISC viscometers to other devices, such as a thermostat, through a PC. Measurement processes and parameters can thus be easily recorded for analysis and documentation. In addition, users can e.g. monitor and control the sample temperature and pH measurement.



designed for scientists

Технические данные

IKA Rotavisc	да
Circulators	да
pH-measurement devices	да
Support of FDA 21 CFR Part 11	да
Labworldsoft Server	да
User management	да
Audit trails	да
Electronic signatures	да
Measurement file database	да
Windows 7	да
Windows 8	да
Windows 10	да

